


Zeiss Axiovert 135 Live Imaging System

Technical Details & Specs



Dr Dale Moulding
UCL Institute of Child Health
Updated Nov 2018

Access to the microscope

- How to book online: <https://ppms.eu/ucl/?ICHFCI>
- How to get training: contact Dale d.moulding@ucl.ac.uk
- Login: Axiovert135 / Password: Axiovert135


Health & Safety

- No food and drink in the room
- CO2 compressed gas cylinder:
 - asphyxiation risk
 - compressed gas risk
- CO2 alarm in the room
- Needles must be kept in the tube provided
- Take your rubbish away with you


The microscope

- Range of objectives available
- Brightfield, phase contrast and fluorescence
- Hamamatsu monochrome camera
- Motorised XY stage and Z focus
- Environmental chamber (temperature and 5% CO₂)

Fluorescence filters




XYZ & shutters/filters control



XY movement

Z Focus



Length gauge to prevent Z drift



Focus Motor on/off

Environmental chamber control

Incubator Temperature controller



Hold to set temperature


CO2 enrichment controller



Environmental chamber



Gas humidifier bottles



- Both bottle caps must closed when in use.
- Open the cap of bottle B when not in use.

Bottle A = filled with clean distilled water up to the 50 mL mark
Bottle B = must be empty

Light sources

Brightfield (Phase contrast)

Green switch back left of scope




Fluorescence

CoolLED pe300 White

Turn on at wall socket.

Display on top of scope will turn on.

All software controlled.



Hamamatsu Monochrome Camera


Hamamatsu Camera controller



Hamamatsu Camera



SPECTRAL RESPONSE CHARACTERISTIC



☒ SPECIFICATIONS

Type number	C4742-80-12AG		
Camera head type	Hermetic vacuum-sealed air-cooled head		
Imaging device	ER-150 progressive scan interline CCD		
Effective number of pixels	1344 (H) ☒ 1024 (V)		
Cell size	6.45 μm (H) ☒ 6.45 μm (V)		
Effective area	8.67 mm (H) ☒ 6.60 mm (V)		
Pixel clock rate	14.75 MHz/pixel		
Frame rate	1 ☒ 1	8.8 frame/s	
	binning	2 ☒ 2	16 frame/s
		4 ☒ 4	27 frame/s
		8 ☒ 8	41 frame/s
Readout noise (r.m.s.) typ.	6 electrons		
Full well capacity typ.	18000 electrons		
Dynamic range* typ.	3000 : 1		
Cooling method	Forced air peltier cooling, with hermetic sealing		
Cooling temperature	- 30 °C		
Dark current	0.03 electrons/pixel/s		
A/D converter	12 bit		
Exposure time	10 μs to 4200 s		
Sub-array	yes		
Contrast enhancement	Analog gain (10times max.) and offset function		
External trigger	yes		
Lens mount	C-mount		
Interface / Output signal (digital output)	IEEE1394-1995 / Non-compressed data (Mono 16)		
External control	IIDC 1394-Based Digital Camera Specification Ver. 1.30		
Line voltage	AC 100 V / AC 117 V / AC 220 V / AC 240 V, 50/60 Hz		
Power consumption	approx. 90VA		
Ambient storage temperature	- 10 °C to + 50 °C		
Ambient operating temperature	0 °C to + 40 °C		
Ambient storage/operating humidity	70 % max. (no condensation)		

Ocular/Camera light path switch

Ocular



Camera



Scale

	Hamamatsu
Axiovert 135	Full Frame bin 1x1
objective	$\mu\text{m}/\text{pixel}$
x10 NA	1.031
x20 NA	0.467
x32 NA	0.318
x40 NA	0.228

Change image unit in Fiji/ImageJ:
Image>Properties ...

Channels (c):

Slices (z):

Frames (t):

Note: c*z*t must equal 1

Unit of length:

Pixel width:

Pixel height:

Voxel depth:

Frame interval:

Origin (pixels):

Global

OK Cancel

Channels (c):

Slices (z):

Frames (t):

Note: c*z*t must equal 1

Unit of length:

Pixel width:

Pixel height:

Voxel depth:

Frame interval:

Origin (pixels):

Global

OK Cancel