


# Start-up of the LSM880 multiphoton confocal



- Make sure an objective is locked into position on the microscope
- Turn on the PC
- Check the key is on (horizontal) - it should always be on (**Picture A**)
- Turn on the MAIN SWITCH (**1**)
- Turn on the SYSTEMS / PC switch (**2**)
- Turn on switch (**3**)



f) Turn on the fluorescence bulb (**Picture B**)

g) Start the Zen software

h) Press "Start System"



i) The lasers are turned on in the software: Go to the Acquisition Tab



j) In the laser window turn on the Argon laser – it takes 5 minutes to warm up. It is outlined in **RED** until it is at full power

k) If using the Multiphoton laser turn this on too (Mai Tai)

# Shutdown of the LSM880 multiphoton confocal

## If there is a user after you:

- Leave everything on
- Transfer your data
- Put the **argon laser on Standby**, turn off the Mai Tai and all other lasers (**Picture C** on previous page)
- Log out of PPMS

## If you are the last user:

- Turn off the lasers in the software (**Picture C**)
- Turn off the fluorescence lamp (**Picture B**)
- Save your data and transfer it using the shared drives
- DO NOT USE USB STICKS OR HARD DRIVES
- Shutdown the software and PC (do not log out of PPMS)
- Turn off switch (3) components (**Picture A**)
- Turn off switch (2) Systems/PC (**Picture A**)
- Wait 5 minutes (until the fan cooling the black box with the laser power supply under the table) has turned off
- After 5 minutes, turn of the Main switch (1) (**Picture A**)
- Check that an objective is in place on the objective turret and half lower it so it is locked in position.