


VENTURA CPAP USER GUIDE

EMERGENCY USE ONLY

**NOT FOR ROUTINE CLINICAL
USE WITHOUT FURTHER
APPROVAL**

**ALL HOSPITALS MUST
WORK WITH THEIR OXYGEN
ENGINEERING TEAMS TO
ASCERTAIN THEIR VIE
OUTFLOW, AND DOWNSTREAM
FLOWS AND PRESSURES TO
SPECIFIC WARD AREAS, BEFORE
DEPLOYING THIS DEVICE**

**ONLY TO BE USED WITH THE
PATIENT CIRCUIT PROVIDED**



USING THE VENTURA CPAP

1. Turn all controls (**A, B, C**) **OFF** (clockwise) and connect the UCL-Ventura to the 4-bar oxygen supply directly or by using an extension hose.
2. Connect the flow generator to the patient circuit shown.
3. Turn the **FLOW ON/OFF** control (**A**) fully **ON** (anticlockwise) and open the **FLOW ADJUSTMENT** control (**B**) (3-4 complete rotations anticlockwise).
4. Increase the oxygen concentration by turning the **OXYGEN ADJUSTMENT** control (**C**) (anticlockwise) to the desired level.
5. Apply mask to the patient with appropriate straps and ensure a good fit.
6. Alter the **FLOW ADJUSTMENT** control (**B**) to ensure optimal flow to:
 - Flutter the mask outlet valve (larger port) throughout the respiratory cycle (the valve will flutter less during inspiration).
 - Hold the blue flow indicator flag (where fitted) horizontal or just below horizontal on inspiration (the flag will remain below horizontal on expiration).
 - A virtually continuous flow at the outlet filter should be felt on the back of a gloved hand throughout the respiratory cycle.
7. Alter the **OXYGEN ADJUSTMENT** control (**C**) if the oxygen concentration has changed.

**If the device becomes damaged it should be
withdrawn from use.**