

## BD FACS Aria III-A Configuration

Laser	Diva Parameter Name	Fluorochromes
488 nm Blue	B 530/30	<b>BB515</b> , FITC, AF488, CFSE, GFP, YFP, TO-PRO-1
	B 586/42	<b>PE</b> , AF555, AF546
	B 616/23	<b>PI</b> , PE/Dazzle 594, PE-CF594, PE-Texas Red, mCherry
	B 695/40	<b>PerCP</b> , PerCP-Cy5.5, PerCP-eFluor 710, PE-Cy5, 7AAD, PI
	B 780/60	<b>PE-Cy7</b>
633nm Red	R 660/20	<b>APC</b> , AF633, AF647, TO-PRO-3
	R 730/45	<b>AF700</b> , APC-Cy5.5, Dye Cycle Red
	R 780/60	<b>APC-Cy7</b> , APC-H7, APC-eFluor780
405nm Violet	V 450/50	<b>BD Horizon 450</b> , Pacific Blue, DAPI, BV421, CFP, PO-PRO-1, Zombie Violet, Cell Trace Violet, BFP
	V 525/50	<b>BD Horizon V500</b> , FVD Aqua Dead, BV510, Pacific Orange
	V 610/20	<b>BV605</b> , Pacific Orange, eFluor605NC
	V 660/20	<b>BV650</b>
	V 710/50	<b>BV711</b> , Pacific Orange, Yellow Dead Cell
	V 780/60	<b>BV785</b>