

Date:	<b>UCL MASS SPECTROMETRY SERVICE</b>	Room No:
		Tel:

Name:	Supervisor: Grant Code:	E-mail:
-------	----------------------------	---------

Sample Filename:_(Starting with your <b>initials then</b> followed by a <b>3</b> figure number):
--

Tick as Required: <input checked="" type="checkbox"/>	EI	CI(pos)	CI(neg)	GC-MS	
No PREFERENCE	APCI(pos)	APCI(neg)	ESI(pos)	ESI(neg)	MALDI-TOF
LC-MS	Accurate Mass	Mass Intensity Listing:		Suitable Solvent:	

Empirical Formula (C <sub>x</sub> H <sub>y</sub> O <sub>z</sub> N):	Mol. Weight (mono-isotopic mass to 5dp):
---	--

Structural Formula:	Toxic: Lachrymator: Smell: Air/Moisture Sensitive:	Please tape sample here
	Impurities Present? YES NO	
	Other Comments:	

Spectra will be sent via e-mail. Please indicate if you want sample returned? YES NO	Email Results
*Please submit sample in solution: MeCN/H <sub>2</sub> O or MeOH/H <sub>2</sub> O at concentration of ~1ng/μL for LC-MS and for ASAP approximate 0.1 mg of powder or a few drops of liquid. Please use vial obtainable from Mass spec.	

Cut here

Date:	<b>UCL MASS SPECTROMETRY SERVICE</b>	Room No:
		Tel:

Name:	Supervisor: Account Code:	E-mail:
-------	------------------------------	---------

Sample Filename:_(Starting with your <b>initials then</b> followed by a <b>3</b> figure number):
--

Tick as Required: <input checked="" type="checkbox"/>	EI:	CI(pos)	CI(neg)	GC-MS	
No PREFERENCE	APCI(pos)	APCI(neg)	ESI(pos)	ESI(neg)	MALDI-TOF
LC-MS	Accurate Mass	Mass Intensity Listing		Suitable Solvent:	

Empirical Formula (C <sub>x</sub> H <sub>y</sub> O <sub>z</sub> N):	Mol. Weight (mono-isotopic mass to 5dp):
---	--

Structural Formula:	Toxic: Lachrymator: Smell: Air/Moisture Sensitive:	Please tape sample here
	Impurities Present? YES NO	
	Other Comments:	

Spectra will be sent via e-mail. Please indicate if you want sample returned? YES NO	Email Results
* Please submit sample in solution: MeCN/H <sub>2</sub> O or MeOH/H <sub>2</sub> O at concentration of ~1ng/μL for LC-MS and for ASAP approximate 0.1 mg of powder or a few drops of liquid. Please use vial obtainable from Mass spec.	