



LSR II Configuration 230718 custom violet

| | | | |
|-------------------|--------|--------------|-----------------------|
| Cytometer: | LSRII | User: | Administrator |
| Cytometer Name: | LSRII | Institution: | N/A |
| Serial Number: | 1 | Software: | BD FACSDiva 8.0.1 |
| Input Device: | Manual | Date: | 9/26/2018 12:26:58 PM |
| Window Extension: | 10.00 | | |




Violet Laser
(405nm)



325 UV Laser
(325nm)

Hoechst Red
Indo-1 (Blue)
BUV[A]496
BUV737
UV[A]-515/30
UV-530/30...




DAPI
Hoechst Blue
Indo-1 (Violet)
BUV[B]496
BUV395
UV[B]-515/30...

Red Laser
(633nm)

APC
AF647
Red-660/20
Red-670/14
To-Pro3

APC-Cy7
APC-H7
Red-780/60



APC-R700
AF700
Red-730/45

Comments: