



LSR II Configuration 260712

Cytometer:	LSRII	User:	Administrator
Cytometer Name:	LSRII	Institution:	ICH
Serial Number:	1	Software:	BD FACSDiva 6.2
Input Device:	Manual	Date:	2/19/2013 1:45:30 PM
Window Extension:	10.00		




Violet Laser
(405nm)



325 UV Laser
(325nm)


Indo-1 (Blue)
UV-530/30
UV-620LP
UV-748ALP



DAPI
Hoechst
Indo-1 (Violet)
UV-405/20
UV-424/44
UV-450/50

Red Laser
(633nm)

APC
Alexa Fluor 647
Red-660/20
To-Pro3



Alexa Fluor 700
Red-730/45

Comments: